

Product Data Sheet

Product Name: DMTr-2'-O-Propargyl-rU-3'-CE-phosphoramidite

Cat. No.: GC81476

Chemical Properties

Cas. No. 171486-62-7


Formula $C_{42}H_{49}N_4O_9P$ M.Wt 784.83

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

DMTr-2'-O-propargyl-rU-3'-CE -Phosphoramidite is a monomeric raw material that can be used for nucleic acid synthesis.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA