

Product Data Sheet

Product Name: MODAG-001

Cat. No.: GC69483

Chemical Properties

Cas. No. 2648041-72-7


Formula $C_{16}H_{15}BrN_4$ M.Wt 343.22

Solubility DMSO : 25 mg/mL (72.84 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA