

Product Data Sheet

Product Name: H-Cit-AMC

Cat. No.: GA22149

Chemical Properties

Cas. No. 93753-78-7

Formula $C_{16}H_{20}N_4O_4$ M.Wt 332.36

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

Cleaved by the cysteine proteases aelurain, a protease from barley related to cathepsin H, and bleomycin hydrolase. See also H-Cit-AMC.HBr.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA