

Product Data Sheet

Product Name: PD-1/PD-L1-IN-63

Cat. No.: GC79374

Chemical Properties

Cas. No. 2691796-92-4


Formula $C_{37}H_{41}Cl_3N_4O_3$ M.Wt 696.11

Solubility DMSO: 104 mg/mL (149.40 mM; Need ultrasonic) Storage Store at 4°C, stored under nitrogen, away from moisture

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

PD-1/PD-L1-IN-63 is an orally active PD-1/PD-L1 inhibitor with an IC_{50} of 9.1 nM. PD-1/PD-L1-IN-63 blocks the PD-1/PD-L1 interaction, induces cancer cell death and inhibits tumor growth. PD-1/PD-L1-IN-63 can be used in the research of melanoma [1].

In Vivo, PD-1/PD-L1-IN-63 (50-100 mg/kg; i.g.; 14 days) exhibits potent dose-dependent antitumor activity in humanized PD-1 C57BL/6J mice inoculated with B16-F10 melanoma [1]. PD-1/PD-L1-IN-63 (50-100 mg/kg; i.g.; 14 days) does not induce cardiotoxicity or myelosuppression in healthy BALB/c

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

mice when administered via oral gavage at doses of 50 or 100 mg/kg once daily for 14 days^[1].

In Vitro, PD-1/PD-L1-IN-63 potently inhibits PD-1/PD-L1 protein interaction with an IC₅₀ of 9.1 nM in an HTRF binding assay^[1]. PD-1/PD-L1-IN-63 binds to both purified human PD-L1 (K_D = 88.3 nM) and murine PD-L1 (K_D = 174 nM) proteins with similar high affinities as measured by ITC^[1]. PD-1/PD-L1-IN-63 (0.05-5 μM) binds dose-dependently to purified human PD-L1 (K_D = 202 nM) and murine PD-L1 (K_D = 471 nM) proteins as measured by SPR^[1]. PD-1/PD-L1-IN-63 (37-1000 nM; 48 h) promotes HepG2 cell death in a dose-dependent manner in a HepG2/hPD-L1 and Jurkat/hPD-1 coculture cell model, increasing HepG2 mortality by 58.2% at 1000 nM^[1]. PD-1/PD-L1-IN-63 (48 h) inhibits the proliferation of MCF-7, H1299, HeLa, HepG2, and Jurkat cells in vitro with IC₅₀ values ranging from 2.3 μM to 4.1 μM^[1].

References:

[1]. Cheng B, et al. Discovery of Novel and Highly Potent Resorcinol Dibenzyl Ether-Based PD-1/PD-L1 Inhibitors with Improved Drug-like and Pharmacokinetic Properties for Cancer Treatment. J Med Chem. 2020 Dec 24;63(24):15946-15959.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA