

Product Data Sheet

Product Name: TLR9-IN-1

Cat. No.: GC70036

Chemical Properties

Cas. No. 2226366-86-3


Formula $C_{23}H_{31}N_7O$ M.Wt 421.54

Solubility DMSO : 500 mg/mL (1186.13 mM; Storage 4°C, protect from light
Need ultrasonic)

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA