

Product Data Sheet

Product Name: Enoxaparin

Cat. No.: GC68034

Chemical Properties

Cas. No. 679809-58-6

Formula

M.Wt

Solubility H₂O : ≥ 100 mg/mL

Storage 4°C, protect from light

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Enoxaparin

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA