

Product Data Sheet

Product Name: Sulfatinib (HMPL-012)

Cat. No.: GC32805

Chemical Properties

Cas. No. 1308672-74-3

SMILES O=S(CC1=CC=CC(NC2=NC=CC(OC3=CC4=C(NC(C)=C4)C=C3)=N2)=C1)(NCCN(C)C)=O


Formula $C_{24}H_{28}N_6O_3S$ M.Wt 480.58

Solubility DMSO : ≥ 30 mg/mL (62.42 mM) Storage Store at $-20^{\circ}C$

General tips For obtaining a higher solubility , please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Protocol

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Kinase experiment: The KDR kinase inhibition activity is tested using the the Z-lyte assay kit. The testing system contains 300 ng/mL of recombinant human KDR catalytic domain, 10 μ M of ATP, 1 μ M of substrate peptide, and a test compound (Sulfatinib) at a series of different concentrations in 384-well plate; total volume is 10 μ L. The enzyme inhibition proceeds at room temperature (25°C), for 1 hour at room temperature on the shaker. 5 μ L of stop solution is added to stop the reaction[2].

Animal experiment: The phamacokinetics of Sulfatinib are studied with male ICR mice (n=6 for each group, weight 20-30g) after a single intravenous and oral dosing at 2.5 and 10mg/kg, respectively. For i.v. dosing formulation, Sulfatinib is dissolved in DMSO (0.25%)-solutol(10%)-ethanol(10%)-physiological saline(79.75%) at the concentration of 0.25 mg/mL. And the p.o. Dosing formulation (1mg/mL) is prepared with 0.5% CMC-Na. After i.v. Or p.o. Dosing, blood samples are collected via the ophthalmic vein at 0 (pre-close), 5, 15, 30 min and 1, 1.5, 2, 4, 8, 24 h, anti-coagulated with heparin-Na. After centrifugation, plasma samples are seprated and protein precipitated with acetonitrilel containing internal standard[2].

References:

[1]. PCT Int. Appl. (2011), WO 2011060746 A1 20110526.

Background

Sulfatinib (HMPL-012) is a potent and highly selective tyrosine kinase inhibitor against VEGFR1/2/3, FGFR1 and CSF1R with IC50s of in a range of 1 to 24 nM.

Sulfatinib inhibits VEGFR1, 2, and 3, FGFR1 and CSF1R kinases with IC50s in a range of 1 to 24 nM, and it strongly blocks VEGF induced VEGFR2 phosphorylation in HEK293KDR cells and CSF1 stimulated CSF1R phosphorylation in RAW264.7 cells with IC50 of 2 and 79 nM, respectively. Sulfatinib also attenuates VEGF or FGF stimulated HUVEC cells proliferation with IC50< 50 nM[1]. Also, it is a hERG inhibitor with IC50 of 6.8 μ M in CHO cell[2].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

In animal studies, a single oral dosing of Sulfatinib inhibits VEGF stimulated VEGFR2 phosphorylation in lung tissues of nude mice in an exposure-dependent manner. Furthermore, elevation of FGF23 levels in plasma 24 hours post dosing suggests suppression of FGFR signaling. Sulfatinib demonstrates potent tumor growth inhibition in multiple human xenograft models and decreases CD31 expression remarkably, suggesting strong inhibition on angiogenesis through VEGFR and FGFR signaling. In a syngeneic murine colon cancer model CT-26, Sulfatinib demonstrates moderate tumor growth inhibition after single agent treatment[1]. After oral dosing of 10 mg/kg, the AUC and Cmax are 397 ng/mL and 138ng/mL in the mouse, respectively[1].

[1]. PCT Int. Appl. (2011), WO 2011060746 A1 20110526.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA