

Product Data Sheet

Product Name: Boc-D-Arg(Tos)-OH

Cat. No.: GA10665

Chemical Properties

Cas. No. 61315-61-5

Chemical Name (2R)-5-[[amino-[(4-methylphenyl)sulfonylamino]methylidene]amino]-2-[(2-methylpropan-2-yl)oxycarbonylamino]pentanoic acid

SMILES CC1=CC=C(C=C1)S(=O)(=O)NC(=NCCCC(C(=O)O)NC(=O)OC(C)(C)C)N


Formula $C_{18}H_{28}N_4O_6S$ M.Wt 428.5

Solubility Soluble in water or 1% acetic acid Storage Desiccate at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

Boc-D-Arg(Tos)-OH

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA