

Product Data Sheet

Product Name: Opemalirsen sodium

Cat. No.: GC79575

Chemical Properties

Cas. No. 2639854-33-2

Formula M.Wt 5480.52

Solubility Storage Store at -20°C, away from moisture

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

Opemalirsen sodium is an antisense oligonucleotide targeted to apolipoprotein L1 (APOL1). It is used for the study of Kidney disorders.

References:

[1]. Mahajan S, Zhao H, Kovacina K, et al. High-sensitivity quantification of antisense oligonucleotides for pharmacokinetic characterization. *Bioanalysis*. 2022;14(9):603-613.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA