

Product Data Sheet

Product Name: Dienestrol

Cat. No.: GC31008

Chemical Properties

Cas. No. 84-17-3

SMILES C/C=C(C1=CC=C(O)C=C1)/C(C2=CC=C(O)C=C2)=C/C


Formula $C_{18}H_{18}O_2$ M.Wt 266.33

Solubility DMSO : 155 mg/mL (581.98 mM) Storage Store at RT

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Dienestrol is a synthetic, non-steroidal estrogen, is an estrogen receptor agonist, for the treatment of menopausal and postmenopausal symptoms.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA