

Product Data Sheet

Product Name: Rusfertide

Cat. No.: GC74315

Chemical Properties

Cas. No. 1628323-80-7


Formula $C_{114}H_{181}N_{27}O_{28}S_2$ M.Wt 2441.95

Solubility DMSO : ≥ 100 mg/mL (40.95 mM) Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37°C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

Rusfertide is a peptide mimetic of natural hepcidin, which targets and degrades ferroportin, reduces serum iron and transferrin-saturation, and thus regulates the production of red blood cells. Rusfertide ameliorates the polycythemia vera, β -thalassemia and hereditary hemochromatosis.

Rusfertide limits the iron toxicity in red blood cells (RBCs) (1 mg/kg, s.c., once every two days, for 49 days) and transferrin-saturation (2.5 mg/kg, s.c., once every two days, for 2 weeks), improves oxygen carrying capacity of RBCs, attenuates the anemia and iron deposition in mice models for β -thalassemia and hereditary hemochromatosis[1].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

References:

[1]. Taranath R, et al., Regulation of iron homeostasis by PTG-300 improves disease parameters in mouse models for beta-thalassemia and hereditary hemochromatosis[J]. Blood, 2019, 134: 3540. [2]. Kremyanskaya M, et al., PTG-300 eliminates the need for therapeutic phlebotomy in both low and high-risk polycythemia vera patients[J]. Blood, 2020, 136: 33-35.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA