

Product Data Sheet

Product Name: Splenopentin

Cat. No.: GA23512

Chemical Properties

Cas. No. 75957-60-7

Formula $C_{31}H_{51}N_9O_9$

M.Wt 693.8

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

Splenopentin, RKEVY, is an immunomodulatory pentapeptide corresponding to the amino acid residues 32-36 of the splenic hormone splenin.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA