

Product Data Sheet

Product Name: GNE-140

Cat. No.: GC79680

Chemical Properties

Cas. No. 1809794-70-4

Formula $C_{25}H_{23}ClN_2O_3S_2$

M.Wt 499.04


Solubility

Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

GNE-140 (Compound 28) is an external racemic mixture of (R)-GNE-140 and (S)-GNE-140. GNE-140 is an effective lactate dehydrogenase (LDH) inhibitor, with the IC_{50} values for LDHA and LDHB being 9 nM and 30 nM respectively. GNE-140 can be used in cancer research [1].

References:

[1]. Purkey HE, et al. Cell Active Hydroxylactam Inhibitors of Human Lactate Dehydrogenase with Oral Bioavailability in Mice. ACS Med Chem Lett. 2016 Aug 26;7(10):896-901.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA