

Product Data Sheet

Product Name: (S)-(+)-3-Aminopiperidine dihydrochloride

Cat. No.: GD00404

Chemical Properties

Cas. No. 334618-07-4

Formula $C_5H_{12}N_2 \cdot 2HCl$

M.Wt 173.08


Solubility

Storage argon filling

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA